

## Urgent Field Safety Notice

**PrisMax** Systems and **TherMax** Blood Warmer Units  
**FA Number:** FAV-2025-005  
**Manufacturer:** Baxter Healthcare SA (CH-MF-000026124)  
**Type of Action:** Correction

01<sup>st</sup> August 2025

Dear Dialysis Provider:

Vantive is issuing a Correction for the **PrisMax** Systems and **TherMax** Blood Warmer Units listed below due to the following three issues:

- **PrisMax** System alarm *T2309: Air Detected in Prime* may be raised after priming with a low liquid level in the deaeration chamber. The system correctly detects the presence of air within the disposable set and raises the alarm at the end of priming if the condition is present. Users should follow the on-screen instructions to reprime the disposable set.
- **PrisMax** System alarm *T0830: Blood Leak Detected* may be raised when no blood leak is present, or the operator may have difficulty normalizing the Blood Leak Detector (BLD) leading to **PrisMax** System alarms *T1313: BLD Normalize Failed*, *T0853: Normalization Failed*, or *T1205: BLD Self-Test Failure*. Operators should follow the on-screen instructions and return the blood in the disposable set if connected to a patient.
- **Thermax** System may be unable to detect the presence of a bag on the **Thermax** Blood Warmer Unit leading to difficulty setting up therapy, or **PrisMax** System alarm *T2284: Thermax Disposable Not Inserted*. Operators should follow the on-screen instructions and return the blood in the disposable set if connected to a patient.

If the operator continues to receive any of these alarms or further assistance is required, please contact your local Vantive Clinical Specialist or local Vantive Account Manager.

Vantive is currently investigating these issues and will be correcting the impacted devices.

### Affected Product

Product Code	Product Description	Serial Numbers
955558	<b>PrisMax</b> V2-ROW	All
955725	<b>PrisMax</b> V3 Control Unit-ROW	All
955515	<b>TherMax</b> Blood Warmer Unit ROW	All

### Hazard Involved

The occurrence of these three issues could lead to a delay or interruption of therapy, and in certain situations could result in blood loss. However, blood return is possible to avoid blood loss. In the event of blood loss, the amount of blood loss would be limited to the amount of blood contained in the disposable set. Vantive is aware that users could decide to not return the blood in accordance with their clinical practices. Patients that are hemodynamically unstable and anemic from the onset are in the high-risk group for blood loss. Vantive has received one complaint of serious injury related to blood loss associated with the above issues.

## Actions to be Taken by Customers

1. Operators may continue to safely use the **PrisMax** Systems and **TherMax** Blood Warmer Units per the associated Operator's Manuals and on-screen instructions. Refer to the enclosed Attachment A for details. If further assistance is required, please contact your local Vantive Clinical Specialist or local Vantive Account Manager.
2. Please ensure that all operators of these devices are made aware of this notification. Vantive recommends that a copy of this notification is posted within the appropriate location of your facility.
3. As corrections become available, a local Vantive representative will contact your facility to determine the correction plan and schedule the associated correction(s) for impacted devices.
4. Complete the enclosed customer reply form and return it to Vantive by scanning and e-mailing it to [Ireland.CORP.FA@vantive.com](mailto:Ireland.CORP.FA@vantive.com), even if you don't have any inventory. Returning the customer reply form promptly will confirm your receipt of this notification and prevent you from receiving repeat notices.
5. If you distributed this product to other facilities or departments within your institution, please forward a copy of this communication to them.

## Further Information and Support

For general questions regarding this communication contact [Ireland.CORP.FA@vantive.com](mailto:Ireland.CORP.FA@vantive.com).

Any adverse reactions or quality problems experienced with the use of these products may be reported using one of the following options:

Reporting product quality complaints:

- Portal: [Vantive - Product Feedback Portal](#)

Reporting adverse events with drugs:

- Email: [UK.Ireland.GPS.AEReporting@vantive.com](mailto:UK.Ireland.GPS.AEReporting@vantive.com)

We apologize for any inconvenience this may cause you and your staff.

Sincerely,



Mike Little  
Head of Sales UKI  
Vantive Health Limited

Enclosures: Vantive Customer Reply Form

Attachments: Attachment A: Operator's Manual and Graphical User Interface (GUI) Guidance

**CUSTOMER REPLY FORM related to Urgent Field Safety Notice (Product Correction)**  
**FAV-2025-005 dated 01st AUGUST 2025**

**Product Name:** PrisMax V2-ROW, PrisMax V3 Control Unit-ROW and TherMax Blood Warmer  
Unit ROW

**Product code:** 955558, 955725 and 955515

**Serial Number: All**

Please complete and return one copy of this form per facility by e-mail ([Ireland.CORP.FA@vantive.com](mailto:Ireland.CORP.FA@vantive.com))  
as confirmation that you have received this notification.

Facility Name and Address:	
Reply Confirmation Completed By ( <i>Please Print</i> ):	
Title ( <i>Please print</i> ):	
Email and/or Telephone Number (including Area Code):	

Your signature below indicates that you have received the attached letter; performed the actions as outlined in the letter as needed; and disseminated this information to staff and other services or facilities as applicable.

<b>Signature/Date:</b>	
REQUIRED FIELD	