



Field Safety Notice
CONMED Corporation
AirSeal® Access Ports

June 24, 2026

Dear Risk Manager/ Safety Coordinator/Purchasing Agent:

CONMED Corporation is sending this communication to provide you with important information concerning the Instructions for Use (IFU) for the AirSeal® Access Ports. This communication affects the catalog numbers listed below for customers who received devices manufactured from 14 June 2024 to 12 Jan 2026. The Lot numbers impacted as a result of this Field Safety Notice is from 202406144 through 202601124.

CONMED is **not** requesting users to return, remove or dispose any of the *Airseal Access Ports* from the market.

Table 1: Catalog number impacted

Catalog No.	Device Identifier	Device Name
iASB5-150	10845854045695	AirSeal® 5 mm Smooth Access Port and Obturator with Blunt Tip, 150mm Length
iASB12-120	10845854045725	AirSeal® 12 mm Access Port and Obturator with Blunt Tip, 120mm Length
iASB12-100	10845854045602	AirSeal® 12mm Access Port and Obturator with Blunt Tip, 100mm Length
iAS8-120LP	10845854045657	AirSeal® 8mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 120mm Length
iAS8-100LP	10845854045718	AirSeal® 8mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 100mm Length
iAS5-75LP	10845854046388	AirSeal® 5mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 75mm Length
iAS5-120LP	10845854045657	AirSeal® 5mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 120mm Length
iAS5-100LP	10845854045626	AirSeal® 5mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 100mm Length
iAS12-150LPI	10845854052006	AirSeal® 12mm Access Port and Palm Grip Obturator with Bladeless Optical Tip, 150mm Length
iAS12-120LPI	10845854045688	AirSeal® 12mm Access Port and Low Profile Obturator with Bladeless Optical Tip, 120mm Length
iAS12-100LPI	10845854045640	AirSeal® 12mm Access Port and Palm Grip Obturator with Bladeless Optical Tip, 100mm Length

Reason for the Notice:

CONMED has received an increase in reports of overpressure alerts with AirSeal Access Ports, which in some cases have resulted in delay of treatment. During a review of complaint data, CONMED identified that a previous manufacturing process update may, under certain conditions, contribute to overpressure alerts. In response, the Instructions for Use (IFU) were clarified to ensure proper device tightening.

In the event of an overpressure condition, AirSeal® is designed to respond automatically via a built-in venting system, which activates when the pressure exceeds the set point by more than 3 mm Hg for more than 3 seconds. When the “overpressure” alert is displayed there is a higher-than-expected pressure within the system between the AirSeal® iFS and the AirSeal® access port but does not represent higher-than-expected pressure within the patient.



Risk to Health:

When the “overpressure” alert is displayed there is a higher-than-expected pressure within the system between the AirSeal® iFS and the AirSeal® access port but does not represent higher-than-expected pressure within the patient. The reported incident may result in the loss of pneumoperitoneum due to the system venting gas to compensate for the overpressure. Loss of pneumoperitoneum may result in loss of field of view, which may lead to iatrogenic injury, more complex or converted procedures, and/or delay of treatment.

Although the erroneous overpressure alerts on the device resulting from the issue described above are not indicative of the actual patient cavity pressure, the overpressure alarm should nevertheless be taken seriously when a true overpressure condition occurs within the patient. *Potential risk associated with a true overpressure could include air embolism, arrhythmias (bradycardia, asystole, or cardiac arrest), pneumothorax, kidney or urinary problems, hypoxia, subcutaneous emphysema, delay to treatment, more complex procedures, and potentially death. CONMED has received no reports of patient or user injury or death associated with these overpressure alerts.*

Actions to be taken by the Users:


A replacement AirSeal disposable should be available prior to the procedure in order to minimize potential delays. Before initiating insufflation, the user should ensure the tubing is properly connected between the AirSeal tube set and Airseal Access Ports. If the connection cannot be securely tightened and the “overpressure” alert cannot be resolved, please replace the affected device.

Detailed instructions on how to connect the AirSeal tube set and Airseal Access Ports are as follows:

Remove the manifold plug from the side of the cannula on the AirSeal Access Port and connect the distal end of the Tri-Lumen Filtered Tube Set with Activated Charcoal Filter to the cannula manifold and tighten.

Note: If the AirSeal Tube Set and Access Port are not appropriately tightened, the pressure reading on the device may not reflect the pressure experienced by the patient. A secure connection can be visually confirmed when one of the white fins on the Tube Set bullseye connector aligns directly to the top of the Access Port (Fig. 15). Looking straight into the bullseye of the Access Port while the Tube Set is connected, one of the fins on the connector piece of the Tube Set should sit directly vertical to the Access Port (Fig. 16).

The figures below were added to the updated Instructions for Use to illustrate proper setup. **The prior Instructions for Use did not include these figures.**



Please review your inventory of AirSeal® Access Ports and place a copy of the enclosed Instructions for Use, Part No. AS-DISP-IFU, Rev. F with any inventory in your possession with an active expiration date.



Please dispose of any other Instructions for Use in your possession. Clarification of the IFU is meant to assist all customers when using the AirSeal® Access Ports correctly.

Copies of the corrected Instructions for Use are also available electronically for each of the AirSeal® Access Ports at eifu.conmed.com/eifu. If you have provided these devices to other facilities or distributors, or you are a distributor, please forward this information as appropriate. It is imperative that all end users of these devices receive this notice.

Please do not return any devices as this is a Field Safety Notice. CONMED is not requesting to remove, dispose or return any of the affected devices. Please complete Appendix A or B and return it with acknowledgement

Please respond to this Field Safety Notice that this letter has been read and understood by using Appendix A for Health Care providers and Appendix B for Importers and Distributors.

The appropriate international competent authorities have been notified of this action.

Follow up:

If you have any questions or requests, please contact the Field Action Support Team at 1-800-448-6506 (8:00am to 7:00pm ET, Monday through Friday) or email accessports2026@conmed.com.

If further information relevant to this issue becomes necessary CONMED will communicate as such. CONMED is dedicated to providing safe and reliable products to our customers and their patients. We are committed to manufacturing products of the highest quality, and we sincerely apologize for any inconvenience this may cause you or your staff.

Sincerely,

Kavita Amin
Regulatory Affairs Manager and Recall Coordinator

Attachments:
AirSeal Access Port IFU: AS-DISP-IFU, Rev. F



**AirSeal® Access Ports
Customer Reply Form**

1. Field Safety Notice (FSN) information	
FSN Reference number*	Pre-filled by manufacturer
FSN Date*	Pre-filled by manufacturer
Product/ Device name*	AirSeal® Access Ports
Product Code(s)	iASB5-150
	iASB12-120
	iASB12-100
	iAS8-120LP
	iAS8-100LP
	iAS5-75LP
	iAS5-120LP
	iAS5-100LP
	iAS12-150LPI
	iAS12-120LPI
iAS12-100LPI	
Batch/Serial Number (s)	N/A

2. Customer Details	
Account Number	
Healthcare Organisation Name*	
Organisation Address*	
Department/Unit	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Customer action undertaken on behalf of Healthcare Organisation	
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.
Print Name*	Customer print name here
Signature*	Customer sign here
Date*	

4. Return acknowledgement to sender	
Email:	accessports2026@conmed.com
Customer Helpline:	1-800-448-6506
Postal Address:	ATTN: Khinthandar Lwin 525 French Road Utica, NY 13502 USA
Deadline for returning the customer reply form*	30 days from receiving the notice

Mandatory fields are marked with *



**AirSeal® Access Ports
Distributor/Importer Reply Form**

1. Field Safety Notice (FSN) information	
FSN Reference number*	Pre-filled by manufacturer
FSN Date*	Pre-filled by manufacturer
Product/ Device name*	AirSeal® Access Ports
Product Code(s)	iASB5-150
	iASB12-120
	iASB12-100
	iAS8-120LP
	iAS8-100LP
	iAS5-75LP
	iAS5-120LP
	iAS5-100LP
	iAS12-150LPI
	iAS12-120LPI
iAS12-100LPI	
Batch/Serial Number (s)	

2. Distributor/Importer Details	
Company Name*	
Account Number	
Address*	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Return acknowledgement to sender	
Email:	accessports2026@conmed.com
Distributor/Importer Helpline:	1-800-448-6506
Postal Address:	ATTN: Khinthandar Lwin 525 French Road Utica, NY 13502 USA
Deadline for returning the distributor/importer reply form*	30 days from receiving the notice

4. Distributors/Importers (Check all that apply)	
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.
<input type="checkbox"/>	I have informed the identified customers of this Field Safety Notice
Print Name*	Distributor/Importer print name here
Signature*	Distributor/Importer sign Here
Date *	

Mandatory fields are marked with *

LOT
202409244
202408284
202411204
202412165
2024121611
2024110810
202511044
202511066
202511055
202508274
202509104
202506194
202508134
202508154
202507294
202505194
202507014
202507024
202503114
202503215
202504014
202412164
202412135
202501036
202501294
202501064
202412065
202411225
202408064
202412145
202501065
202502275
202509244
202507254
202508054
202508194
202505064
202504294
202504155
202504304
202501085
202501146
202501274
202501224
202410314
202410315
202407015
202510144
202510094
202506244
202509234
202504074
202412175
202408235
202407294
202511075
202511244
202509154
202509264
202508084
2025080611
202507284
202507264
202504174
202506274
2025040310
202504244

202505144
202501094
202504164
202502214
2024112110
202410245
202410285
202411044
2024112011
202412105
202411266
202512164
202512174
202509174
202510284
202509304
202508074
202508184
202507304
202506204
2025062010
202506174
202506054
2025060510
202507174
202507184
202503296
202503314
2025042910
202505074
2025061010
202506074
202412026
202411306
202411284
202503035
202502245
202412044
202502284
202410144
202410114
202408306
202411114
2024121710
2024121711
202501035
202412195
202412054
202510014
202407014
202509034
202510274
202510184
202406226
202406216
202407024
202407136
202408066
2024071212
202511184
202511054
202406275
202406146
202406175
202409125
202409135
202407184
202407185

202407195
202407194
2025100810
2025100710
202511204
2025092610
202506124
202411285
202507114
202505234
202501164
202410304
202410305
202409065
2025072910
2025072310
202504114
2025050810
202501075
202411134
202411145
202407025
202409024
2025121810
202506094
202411074
2025111110
202508224
2025071410
202507144
202503284
2025060310
202505104
2025050910
202504284
202412035
202501176
202502265
202501046
202502125
202411156
2024121610
2024111111
2024121811
202411185
202410165
202409305
202410105
202410024
2024110710
202412185
202411085
202412055
202407036
202407265
202409205
202407165
202409124
202409025
202407094
2025042210
202505084
202512024
202510024
202512224
2025100610
202505274
202410155

202406205
202508234
202511124
202509294
2025070710
202503194
202505094
202502054
202503034
2024111510
202410154
202512054
202407224
202506064
202512064
202512194
202510304
202511154
2025102310
2025052710
202502285
2024121410
202501175
202410175
202410084
202409175
202406156
202503276
202512094
202509084
202411065
202511064
202511065
202509274
202508214
202510254
202503204
202504044
202503075
202501244
202501154
202407026
202408194
202407116
202512034
202510214
202511014
202510314
202411115
202407084
202507224
202507244
202508124
202411286
202512084
202511194
202503265
202412196
202502145
202409204
202506104
202501024
202511214
2025071010
202505164
202505204
202506254
202503144

202409106
202511114
202508114
2025070210
202506184
2025060410
202506034
202506044
202503055
202501306
202411064
202410094
202407045
202406196
202407315
202407066
202409114
202508044
202510064
202509184
202510074
2025060910
202505096
202501225
202510224
2025082110
202504116
202502016
202410025
202502274
202510294
202511144
202407254
202412174
202507154
202502195
202502045
202411254
202409104
202409245
202406214
2025101510
202408224
202406215
202501166
202512204
202511254
202503115
202510134
202503145
202412184
202407096
202406264
202409174
202512124
202408304
202504144
202502254
202510154
202412186
202505124
202408215
202510244
202511104
202508144
202510104
202507034
202505054

202412064
202408025
202408165
202410225
202408135
202412146
202511074
202509124
2025100310
2025041710
202501084
202406255
202508254
202408085
202407126
202408195
202504095
2024121011
202409194
202409185
202409184
202407054
202502114
202510034
202412194
202507084
202507074
202501135
202412205
202412204
202408084
202411295
202412106
202410254
202408234
2025082210
2025072510
2025072410
202503045
202505284
202501295
202411264
202503184
202410085
202509054
2025102210
202504084
202503196
202504145
202408164
202411055
202411214
202511264
202512184
202509044
202507044
202507045
2025060610
202503255
202503254
202503185
202505224
202505134
2025052110
202505214
202501314
202501216
202501186

202408276
202408226
202411294
2024120411
202412036
202407276
202407305
202407306
202407314
202407044
202407046
202502255
202408134
202507234
202406195
202411124
202508204
202510084
202511134
202411275
2025052210
202503125
202504154
202507164
2025042810
2025040410
202503275
202501174
202407034
202409304
202507314
202502174
202501104
202501074
202411054
202509194
202408014
202410295
2025102410
202510164
202501205
202408094
202505264
2025081110
202407155
202412134
202509114
202503064
202406256
2025121710
202512104
202512154
202509094
202509224
202511046
202507214
2025072810
2025080610
202508064
202506234
2025062310
202503105
202503286
202504254
2024111511
202502184
202502134
202502175

202502135
2024120410
202411084
202411075
202410035
2024121810
202411144
202407035
202409036
202407235
2024072510
202407245
202407255
202407095
202407164
202504034
202501165
202408225
202507104
202506164
202412124
202501144
2025082010
202509254
202504046
2025040411
202501235
202501256
2025092910
202501296
202406276