

Summary of Product Characteristics

1 NAME OF THE MEDICINAL PRODUCT

Nexium 40 mg gastro-resistant tablets

2 QUALITATIVE AND QUANTITATIVE COMPOSITION

Each gastro-resistant tablet contains 44.5 mg esomeprazole as magnesium trihydrate equivalent to 40 mg esomeprazole.

Excipient(s) with known effect

Each gastro-resistant tablet contains 30mg sucrose.

For the full list of excipients, see section 6.1.

3 PHARMACEUTICAL FORM

Gastro-resistant tablet

Product imported from Greece

A pink, oblong, biconvex tablet with '40mg' on one side and  on the other.

4 CLINICAL PARTICULARS

As per PA2242/013/003

5 PHARMACOLOGICAL PROPERTIES

As per PA2242/013/003

6 PHARMACEUTICAL PARTICULARS

6.1 List of excipients

Glycerol monostearate 40-55
Hyprolose
Hypromellose
Iron oxide (reddish-brown) (E 172)
Magnesium stearate
Methacrylic acid ethyl acrylate copolymer (1:1) dispersion 30 per cent
Cellulose microcrystalline
Synthetic paraffin
Macrogols
Polysorbate 80
Crospovidone
Sodium stearyl fumarate
Sugar spheres (sucrose and maize starch)
Talc
Titanium dioxide (E 171)
Triethyl citrate

6.2 Incompatibilities

Not applicable.

6.3 Shelf life

The shelf life expiry date for this product shall be the date shown on the packaging of the product on the market in the country of origin.

6.4 Special precautions for storage

Do not store above 30 °C.

Store in the original package (blister) in order to protect from moisture.

6.5 Nature and contents of container

Blister packs in wallet and/or carton of 28 or 30 Tablets in an over labelled outer carton.

Not all pack sizes may be marketed.

6.6 Special precautions for disposal

No special requirements for disposal.

Administration through gastric tube:

1. Put the tablet into an appropriate syringe and fill the syringe with approximately 25 ml water and approximately 5 ml air. For some tubes, dispersion in 50 ml water is needed to prevent the pellets from clogging the tube.
2. Immediately shake the syringe for approximately 2 minutes to disperse the tablet.
3. Hold the syringe with the tip up and check that the tip has not clogged.
4. Attach the syringe to the tube whilst maintaining the above position.
5. Shake the syringe and position it with the tip pointing down. Immediately inject 5-10 ml into the tube. Invert the syringe after injection and shake (the syringe must be held with the tip pointing up to avoid clogging of the tip).
6. Turn the syringe with the tip down and immediately inject another 5-10 ml into the tube. Repeat this procedure until the syringe is empty.
7. Fill the syringe with 25 ml of water and 5 ml of air and repeat step 5 if necessary to wash down any sediment left in the syringe. For some tubes, 50 ml water is needed.

7 PARALLEL PRODUCT AUTHORISATION HOLDER

Merit Pharmaceuticals Limited
Unit C4/C3
Metropoint Business Park
Kettles Lane
Swords
Co Dublin
K67 RH92
Ireland

8 PARALLEL PRODUCT AUTHORISATION NUMBER

PPA23080/003/002

9 DATE OF FIRST AUTHORISATION/RENEWAL OF THE AUTHORISATION

Date of first authorisation: 28th May 2010

10 DATE OF REVISION OF THE TEXT

February 2025