

# Summary of Product Characteristics

## 1 NAME OF THE MEDICINAL PRODUCT

Prolastin 5000 mg, powder and solvent for solution for infusion

## 2 QUALITATIVE AND QUANTITATIVE COMPOSITION

One vial of powder contains approx. 5000 mg alpha<sub>1</sub>-proteinase inhibitor, human\*.

After reconstitution with 200 mL solvent, the solution contains approx. 25 mg/mL alpha<sub>1</sub>-proteinase inhibitor (human).

\*Produced from the plasma of human donors

### Excipients with known effect:

Prolastin contains 2.76 mg sodium per mL of reconstituted solution (120 mmol/L).

For the full list of excipients, see section 6.1.

## 3 PHARMACEUTICAL FORM

Powder and solvent for solution for infusion

Powder or friable mass: White or pale yellow or pale brown.

Solvent: Clear colourless solution.

The reconstituted solution appears clear to slightly opalescent, colourless, pale green, pale yellow or pale brown.

## 4 CLINICAL PARTICULARS

### 4.1 Therapeutic indications

Prolastin is indicated for long-term augmentation therapy in subjects with documented severe alpha<sub>1</sub>-proteinase inhibitor deficiency (e.g. genotypes PiZZ, PiZ(null), Pi (null,null) and PiSZ). Patients are to be under optimal pharmacologic and non-pharmacologic treatment and show evidence of progressive lung disease (e.g. lower forced expiratory volume per second (FEV1) predicted, impaired walking capacity or increased number of exacerbations) as evaluated by a healthcare professional experienced in the treatment of alpha<sub>1</sub>-proteinase inhibitor deficiency.

### 4.2 Posology and method of administration

Physicians experienced with chronic obstructive lung diseases should initiate the treatment and supervise the first infusions. Subsequent infusions can be administered by a healthcare professional, see section 4.4.

The duration of treatment is at the discretion of the attending physician. A specific limit on the duration of treatment is not envisaged.

### Posology

#### *Adults, including elderly patients*

Unless otherwise prescribed, a once-weekly dose of 60 mg active ingredient/kg body weight (equivalent to 180 ml reconstituted solution for infusion containing 25 mg/ml alpha<sub>1</sub>-proteinase inhibitor (human) in the case of a patient weighing 75 kg) as a short-term infusion is usually sufficient to keep the serum alpha<sub>1</sub>-proteinase inhibitor level constantly over 80 mg/dl which correlates with pulmonary levels of 1.3 microM. Theoretically this serum and epithelial-lining fluid levels are estimated protective levels against further worsening of the pulmonary emphysema.

#### *Paediatric population*

The safety and efficacy of Prolastin in the paediatric population (<18 years) has not yet been established. No data are available.

### Method of administration

Prolastin should only be administered by intravenous infusion after reconstitution.

The dry powder must be dissolved with the co-packed solvent (water for injections) as described in section 6.6 and administered using a suitable infusion set (not included).

The prepared solution must be used within 3 hours of its preparation.

The infusion rate should not exceed 0.08 ml/kg body weight per minute (equivalent to 6 ml per minute in a patient weighing 75 kg). This infusion rate may be adjusted based upon patient tolerability.

For instructions on reconstitution of the medicinal product before administration, see section 6.6.

### 4.3 Contraindications

Prolastin must not be used in patients with

- selective IgA deficiency who are known to have antibodies against IgA, as allergic reactions to the point of anaphylactic shock may occur in such cases,
- known hypersensitivity to alpha<sub>1</sub>-proteinase inhibitors or to any of the excipients listed in section 6.1 (see also section 4.4).

### 4.4 Special warnings and precautions for use

The recommended infusion rate given under section 4.2 should be adhered. If any reaction takes place that might be related to the administration of Prolastin, the rate of infusion should be decreased, or the administration should be stopped, as required by the clinical condition of the patient.

Since Prolastin can cause a transient increase in blood volume, particular caution is necessary in patients with severe heart failure and patients at risk of circulatory overload.

#### Hypersensitivity

Rarely hypersensitivity reactions may occur, even in patients who have tolerated previous treatment with human alpha<sub>1</sub>-proteinase inhibitor. In the event of a severe hypersensitivity reaction (with a fall in blood pressure to < 90 mmHg, dyspnoea or even anaphylactic shock), Prolastin should be discontinued immediately and suitable therapy, with treatment for shock as necessary, should be instituted.

#### Home-treatment

There are limited data regarding the use of Prolastin in home-treatment.

Potential risks associated with home-treatment are related to the handling and administration of the medicinal product as well as to the handling of adverse reactions. Patients should in any case be informed of signs of hypersensitivity reactions.

The decision of whether a patient is suitable for home-treatment is made by the attending physician, who should ensure appropriate training is provided (e.g. regarding reconstitution, use of transfer device for reconstitution, assembly of intravenous tubing, infusion techniques, maintenance of a treatment diary, identification of adverse reactions and measures to be taken in case such reactions occur) and the use is reviewed at regular intervals.

#### Transmissible agents

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens.

The measures taken are considered effective for enveloped viruses such as human immunodeficiency virus (HIV), hepatitis B virus (HBV) and hepatitis C virus (HCV). The measures taken may be of limited value against non-enveloped viruses such as hepatitis A and parvovirus B19.

Parvovirus B19 infection may be serious for pregnant women (fetal infection) and for individuals with immunodeficiency or increased erythropoiesis (e.g. haemolytic anaemia).

Appropriate vaccination (hepatitis A and B) should be considered for patients in regular repeated receipt of human plasma-derived proteinase inhibitors.

Traceability

Every time that Prolastin is administered to the patient, the name and batch number should be clearly recorded in order to maintain a link between the patient and the batch of the product.

Smoking

Prolastin therapy cannot be denied to smokers. However, since the efficacy of Prolastin will be compromised by tobacco smoke in the lungs, cessation of smoking is strongly recommended.

Sodium content

This medicinal product contains approximately 552.0 mg (24 mmol) sodium per 5000 mg vial. In case of a patient of 75 kg body weight the sodium content of the recommended dose is equivalent to 24.84 % of the WHO recommended maximum daily intake of 2 g sodium for an adult.

This should be taken into account for patients on a controlled sodium diet.

**4.5 Interaction with other medicinal products and other forms of interaction**

No interaction studies have been performed.

**4.6 Fertility, pregnancy and lactation**Pregnancy

For Prolastin no clinical data on exposed pregnancies are available. Animal studies have not been conducted. Caution should be exercised when prescribing to pregnant women.

Breast-feeding

It is unknown whether alpha<sub>1</sub>-proteinase inhibitor is excreted in human breast milk. The excretion of alpha<sub>1</sub>-proteinase inhibitor in milk has not been studied in animals. A decision on whether to continue/discontinue breast-feeding or to continue/discontinue therapy with Prolastin should be made taking into account the benefit of breast-feeding to the child and the benefit of Prolastin therapy to the woman.

Fertility

No animal fertility studies have been conducted with Prolastin. Since human alpha<sub>1</sub>-proteinase inhibitor is an endogenous human protein; thus, no adverse effects on fertility are expected when given at recommended doses.

**4.7 Effects on ability to drive and use machines**

Prolastin has no or negligible influence on the ability to drive and use machines.

**4.8 Undesirable effects**Summary of the safety profile

Treatment with Prolastin may result in known reactions as fever, flu-like symptoms, dyspnoea, urticaria, nausea, etc. Uncommon or rare immunological reactions may occur as with any protein treatment, even when the patient has shown no hypersensitivity or allergic reaction to previous administration. This would include allergic reactions such as urticaria or dyspnoea, very rarely anaphylaxis (see section 4.4). Symptoms which are of possible immunological origin should be evaluated before patients are re-challenged with therapy.

Tabulated list of adverse reactions

The table presented below is according to the MedDRA system organ classification (SOC and Preferred Term Level). Frequencies have been evaluated according to the following convention: very common ( $\geq 1/10$ ); common ( $\geq 1/100$  to  $< 1/10$ ); uncommon ( $\geq 1/1,000$  to  $< 1/100$ ); rare ( $\geq 1/10,000$  to  $< 1/1,000$ ); very rare ( $< 1/10,000$ ), not known (cannot be estimated from the available data).

The following undesirable effects have been observed during treatment with Prolastin:

System organ class	Uncommon	Rare	Very rare
--------------------	----------	------	-----------

	> 0.1% to <1%	>0.01% to <0.1%	< 0.01%
<b>Immune system disorders</b>	Urticaria	Hypersensitivity reactions	Anaphylactic shock
<b>Nervous system disorders</b>	Dizziness / dazed state headache		
<b>Cardiac disorders</b>		Tachycardia	
<b>Vascular disorders</b>		Hypotension Hypertension	
<b>Respiratory, thoracic and mediastinal disorders</b>	Dyspnoea		
<b>Skin and subcutaneous tissue disorders</b>	Rash		
<b>Gastrointestinal disorders</b>	Nausea		
<b>Musculoskeletal and connective tissue disorders</b>	Joint pain / arthralgia	Back pain	
<b>General disorders and administration site conditions</b>	Chills, fever, flu-like symptoms, chest pain		

For information on viral safety, see section 4.4.

#### Reporting of suspected adverse reactions

Reporting suspected adverse reactions after authorisation of the medicinal product is important. It allows continued monitoring of the benefit/risk balance of the medicinal product. Healthcare professionals are asked to report any suspected adverse reactions via HPRA Pharmacovigilance, Website: [www.hpra.ie](http://www.hpra.ie)

#### **4.9 Overdose**

Consequences of overdose are not known.

In the event of overdose, the patient should be observed closely for the occurrence of undesirable effects and supportive measures should be available as necessary.

### **5 PHARMACOLOGICAL PROPERTIES**

#### **5.1 Pharmacodynamic properties**

Pharmacotherapeutic group: proteinase inhibitor, ATC code: B02AB02

Alpha<sub>1</sub>-proteinase inhibitor is a normal constituent of human blood that inhibits the activity of neutrophil elastase, amongst other enzymes. Alpha<sub>1</sub>-proteinase inhibitor has a molecular weight of 51kDa and belongs to the family of serine protease inhibitors.

#### Mechanism of action

It is currently assumed that the pathogenesis of emphysema in alpha<sub>1</sub>-proteinase inhibitor deficiency is attributable to a chronic biochemical imbalance between elastase and alpha<sub>1</sub>-proteinase inhibitor. Elastase, which is synthesised by pro-inflammatory cells in the lower respiratory tract, is capable of breaking down elastic tissue. One of the principal inhibitors of elastase is alpha<sub>1</sub>-proteinase inhibitor, which is lacking in hereditary alpha<sub>1</sub>-proteinase inhibitor deficiency. As a result, the alveolar structures remain unprotected against the elastase that is released by the neutrophils in the lower respiratory tract and to which they are therefore chronically exposed.

This leads to progressive degradation of the elastic tissue and when serum alpha<sub>1</sub> antitrypsin levels decrease below 80 mg/dl, (which correlates with pulmonary levels of 1.3 microM), this is associated with an increased risk for development of emphysema.

#### Pharmacodynamic effects

The administration of Prolastin increases and maintains serum levels and lung epithelial lining fluid (ELF) levels of alpha<sub>1</sub>-proteinase inhibitor leading to a slowdown of the progression of emphysema.

In two controlled observational registries, the most significant slowing of reduction in the rate of FEV1 has been observed in patients with FEV1 35 to 60 % of predicted.

#### **5.2 Pharmacokinetic properties**

After intravenous administration, virtually 100 % of the alpha<sub>1</sub>-proteinase inhibitor dose is immediately available in the patient's bloodstream. The mean in vivo recovery rate is 4.2 mg/dl per kg body weight. The in vivo half-life is approximately 4.5 days.

Pharmacokinetic parameters evaluated for Prolastin in the ChAMP study, a multi-center, randomized, double-blind, crossover trial to evaluate the pharmacokinetic comparability of Alpha-1 MP(Modified Process) to Prolastin in 24 adult subjects with Alpha1-antitrypsin (AAT) deficiency are presented in Table 1 below:

Table 1: Primary PK Endpoint and Other Key PK Parameters Determined from the Potency Assay Results

Treatment	AUC <sub>0-7 days</sub> Mean (%CV)	C <sub>max</sub> Mean (%CV)	t <sub>max</sub> Median (Range)	t <sub>1/2</sub> Mean (%CV)	Mean C <sub>trough</sub> Mean (%CV)
	(mg*h/ml)	(mg/ml)	(hr)	(hr)	(mg/ml)
Prolastin	152.4 (16%)	1.848 (15%)	0.820 (0.25-2.90)	139.3 (18%)	0.574 (20%)

### 5.3 Preclinical safety data

The active ingredient of Prolastin, alpha<sub>1</sub>-proteinase inhibitor, is obtained from human plasma and behaves like endogenous plasma constituents. Administration of a single dose of Prolastin to various animal species, as well as administration of daily doses during five consecutive days to rabbits, showed no indications of toxic effects. Additional preclinical studies with repeated dosing (chronic toxicity, carcinogenicity, reproduction toxicity) were not conducted. These studies cannot usefully be performed in conventional animal models, as antibodies are expected to be formed as a result of administration of heterologous human proteins. Since human alpha1-proteinase inhibitor is a protein and a physiological constituent of human blood, it is not expected to present carcinogenic, genotoxic, or teratogenic effects.

## 6 PHARMACEUTICAL PARTICULARS

### 6.1 List of excipients

Powder: Sodium chloride, Sodium dihydrogen phosphate  
Solvent: Water for injections

### 6.2 Incompatibilities

In the absence of compatibility studies Prolastin must not be mixed with medicinal products or other solutions for infusion.

### 6.3 Shelf life

3 years

The reconstituted solution should always be used within 3 hours of its preparation.

Once prepared, the solution for infusion should not be stored in a refrigerator. Discard any unused solution according to local requirements.

### 6.4 Special precautions for storage

Do not store above 25°C.

Do not freeze.

For storage conditions after reconstitution of the medicinal product, see section 6.3.

### 6.5 Nature and contents of container

Powder:

Type II glass vial with chlorobutyl rubber stopper and aluminium cap.

Solvent:

Type II glass vial with bromobutyl rubber stopper and aluminium cap.

### **Pack sizes**

Single pack of Prolastin 5000 mg, powder and solvent for solution for infusion contains:

One powder vial (5000 mg alpha<sub>1</sub>-proteinase inhibitor, human),

one solvent vial (200 mL water for injections),

one transfer device for reconstitution.

### **6.6 Special precautions for disposal and other handling**

The dry powder is to be brought into contact with and dissolved in the contents of the co-packed vial containing water for injections, as described below. The reconstituted solution appears clear to slightly opalescent, colourless, pale green, pale yellow or pale brown. Total reconstitution should be obtained within 15 minutes.

### **Preparation of the reconstituted solution for infusion**

1. Use aseptic technique (clean and sanitized) in order to maintain sterility. Perform the reconstitution procedure on a flat work surface.
2. Ensure that the vials of Prolastin powder and the solvent (Sterile Water for Injections) are at room temperature (20-25°C).
3. Remove the protective cap from both the Prolastin vial and the solvent vial and clean the collar rims and the stoppers with an alcohol swab. Allow the rubber stoppers to dry.

4. Open the sterile transfer device package by peeling away the lid completely. Do not remove the device from the package.




5. Place the solvent vial upright on the even surface and hold it securely. Without removing the outer package press the blue end of the transfer adapter straight down until the spike penetrates the stopper and snaps. Avoid rotating.



6. Remove the clear outer packaging from the transfer adapter and discard it.



7. Place the Prolastin powder vial upright on the surface. Turn the unit consisting of the adapter and the solvent vial upside down by 180°. Push it with the clear/white end of the adapter straight down – without rotating – until the spike penetrates the stopper and snaps.



8. Due to the vacuum in the powder vial the solvent transfer will start automatically. Wait for complete transfer of the solvent. Remove the adapter with the connected solvent vial in an approximate angle of 45°.



9. Gently swirl the Prolastin vial until the powder is completely dissolved. Do not shake to avoid foaming. Do not touch the stopper. Proceed with the administration of the product by aseptic technique.



10. If more than one vial of product will be needed to achieve the required dose, repeat instructions above using the additional package containing a new transfer adapter. Do not reuse the adapter.

Only clear to slightly opalescent, colourless, pale green, pale yellow or pale brown solutions, free from visible particles should be used. The reconstituted solution must always be used within 3 hours of its preparation. Any unused medicinal product or waste material should be disposed of in accordance with legal requirements.

## 7 MARKETING AUTHORISATION HOLDER

Grifols Deutschland GmbH  
Colmarer Strasse 22  
60528 Frankfurt  
Germany

## 8 MARKETING AUTHORISATION NUMBER

PA1405/002/003

## 9 DATE OF FIRST AUTHORISATION/RENEWAL OF THE AUTHORISATION

Date of first authorisation: 5th January 2007  
Date of latest renewal: 21st March 2031

## 10 DATE OF REVISION OF THE TEXT

January 2026